



# KaVo QUATTROcare®

Automatic Handpiece Maintenance and Purging System  
- Ideal handpiece care... simplified!

**100%  
KaVo.**  
First Choice, Best Value!

**2004**

**REALITY**  
Five Star Award




# HandPeace

The fast and convenient solution  
to protect your handpiece  
investment!

The correct maintenance of dental handpieces is critical to ensure top performance and longevity. Proper care means reduced repair costs, less frustration and the elimination of costly downtime.

Proper cleaning of rotating instruments is essential before lubrication and sterilization. If debris remains in the inside of the instrument, it can cause handpiece failure.

## Don't lubricate the debris in your handpieces!




The KaVo QUATTROcare ensures proper cleaning and removal of debris. The unique cleaning and maintenance mixture expands 200-300 times inside the instrument producing a highly effective cleaning foam. The automatic dispenser delivers the precise amount of cleaner and lubricant time after time leading to longer handpiece life and less repair needs. QUATTROcare has proven itself in thousands of dental offices all over the world and is rated #1 by REALITY and The Dental Advisor. The QUATTROcare also received a 2003 Townie Choice Award.

*"Simply put, the QUATTROcare is the best thing I ever did for my handpieces!!! We had been routinely repairing handpieces every 3 or 4 months, but after purchasing the QUATTROcare we have not needed one repair - that's none, zero, zippo!! The entire staff is absolutely thrilled with it. No more guess work when lubricating or purging. What's not to love?! It's a simple, fast, clean machine that does more than it promises..."*

*-Michael A. Esposito, DDS, Stroudsburg, Pennsylvania*



Here's how simple it is:



Position the instruments...



...close the door...



...press Start...



...and everything's done...

# KaVo QUATTROcare

## Load.

## Start.

## That's it!

KaVo QUATTROcare offers the following advantages and benefits:

### EASY:

- Designed to clean, lubricate and purge your instruments with the push of a button - **simplifying maintenance**
- Automatic care and maintenance procedure
- Easy loading and operation - **user friendly**

### FAST:

- **Up to four (4) handpieces** or air motors are cleaned, lubricated and purged in the same cycle
- Care and maintenance of one complete set of handpieces in the same cycle is possible - **no need to change adapters** between high speeds and low speed attachments
- Assistant **saves valuable time**, can concentrate on other aspects of work
- Flexible programming allows you to choose and set total cycle time from 20-45 seconds per handpiece
- Handpieces can be sterilized immediately after maintenance cycle - no need for separate purging of the handpiece
- Leads to **improved organization** in the dental office


### ECONOMICAL:

- QUATTROcare **increases the longevity** of your handpieces
  - Less wear of handpieces, less repairs
  - Smoother running of handpieces
  - Protects your handpiece investment
- QUATTROcare **provides an excellent ROI (Return On Investment)** due to saving valuable staff time and less repairs
- QUATTROcare can be used for virtually all manufacturers' high speeds, low speed contra-angles, straight handpieces, air driven scalers and air motors


### SAFE:

- Front-loading enclosed maintenance chamber - hygienically isolates entire process
- QUATTROcare **eliminates mishandling or insufficient care** and maintenance
- Completely controlled cleaning and purging procedure
- Visual as well as audio signal for error indication
- Indicator shows when to replenish spray

# KaVo QUATTROcare specifications:



**KaVo QUATTROcare**  
itself is quickly and  
easily cleaned.



**KaVo QUATTROcare...**  
**LOAD. START. THAT'S IT.**

**KaVo QUATTROcare is easy to use –**  
handpiece maintenance is so simple with  
QUATTROcare!

## KaVo QUATTROcare Unit and Accessories

Part	Order No.
KaVo QUATTROcare Unit	0411.7800
KaVo MULTIflex adapter (for Cellular Optic High Speeds, SONICflex airscaler, INTRAFlex low speeds)	0411.7991
KaVo E-type adapter (for INTRAMatic low speed and surgical attachments)	0411.7981
KaVo INTRA-head adapter	0411.7941
Midwest "Shorty" attachment adapter	1000.8459
Bien Air High Speed adapter	1001.7860
W&H High Speed adapter	1001.7888
NSK High Speed adapter	1000.8766
KaVo QUATTROcare spray (6 cans)	0411.7640

## Spacesaving countertop design –

QUATTROcare fits virtually anywhere.

### QUATTROcare dimensions:

Height: 15"  
Width: 13-9/16"  
Depth: 11"

Width with door opened (right hand side –  
to exchange spray can): 16-15/16"  
Depth with flap open (front): 18-7/8"

Weight: 22 lbs.

### QUATTROcare requirements:

- 120V, 50/60Hz, 20 amps electrical connection
- Simply connects to 1/4" airline  
(58 to 145 psi air pressure – built in  
automatic pressure regulator controls flow)

**KaVo America Corporation**  
340 East Main Street, Lake Zurich, IL 60047  
1.800.KaVo USA (1.800.323.8029) • Fax 1.847.550.6825  
www.kavousa.com



KaVo. Dental Excellence.