

NATURAL ELEGANCE® NANO-HYBRID COMPOSITE

**SAVE TIME AND
DELIVER LONG-
LASTING, NATURAL
LOOKING RESULTS**



For more information, see product Instructions for Use.

Order today at [HENRYSCHEIN.NL](https://www.henryschein.nl)



EXCLUSIVE BRANDS
& PRODUCTS



QUALITY YOU
CAN RELY ON



SMART
SAVINGS

NATURAL ELEGANCE® NANO-HYBRID COMPOSITE

NATURAL Elegance NH is a nano-hybrid composite for everyday dentistry. This universal composite is completely focused on simplifying shade selection without neglecting the aesthetic aspects. One duo shade matches 2 Vita™ shades (e.g., A1/B1 or A2/B2). As a result, a small number of shades is sufficient to cover a wide range of tooth shades.



Advantages

Firm, non-sticky thixotropic formula for ideal sculpting and control

Excellent margin adaptation adheres to cavity wall without leaving micro-holes

Duo shade system reduces the number of colors required and helps achieve natural results

Long-lasting wear resistance

Highly polishable HCM (high cohesive molecular)

Radiopaque

Universal composite appropriate for anterior and posterior restorations of all classes

Excellent optical self-blending properties

NATURAL ELEGANCE® NANO-HYBRID COMPOSITE

Syringe, 4g

A1/D2	9003538
A2/B2	9003539
A3/D3	9003540
A3,5/B3	9003541

Tips, 20 x 0.25g

A1/D2	9003545
A2/B2	9003546
A3/D3	9003547
A3,5/B3	9003548

For more information, see product Instructions for Use.

Order today at **HENRYSCHEIN.NL**



EXCLUSIVE BRANDS
& PRODUCTS



QUALITY YOU
CAN RELY ON



SMART
SAVINGS