



EDGEAPEX HP™ EDGEFLOW™ EDGEPACK™ EDGEAPEX™

# EDGEPOWER

MORE POWER FOR ENDODONTICS

The Edge Power portfolio was designed to work seamlessly with our **NiTi and Obturation solutions.**





# EDGEAPEX HP™

**EDGEAPEX HP™** Designed to enhance precision and efficiency in endodontic procedures. Cordless and with an Integrated Apex Locator for maximum convenience.

## Benefits

- Compatible with the majority of the existing file systems for better adaptability.
- Exceptional versatility with two movement options, continuous rotation and reciprocation.
- Cordless portable motor with wireless charging, auto on/off, and an integrated Apex Locator for maximum convenience.
- High-performance brushless motor with quiet operation and an extended service life.
- It adopts real-time feedback technology and dynamic torque contra angle, effectively preventing file separation.

### Specifications:

- Great flexibility with a 360° rotation contra angle.
- Ultra-small head and slim neck for enhanced visibility and improved access to difficult areas.
- 10 available memories and pre-program settings for EdgeEndo file systems to save time.
- 4 available motions: CW / CCW / REC / ATR
- Multiple Apical functions
- Right-handed and left-handed adjustable display
- RPM range: 100 - 2500
- Torque range: 0.4 - 5.0 Ncm



Item Code: **EAPEXHP**  
 Motor Hand piece: **EAPEXRU**  
 Charging Base: **EAPEXCB**

## Accessories included in package

The box includes a motor handpiece, a charging base, and a contra angle. Measuring Wire, Lip Hook (x2), File Clip (x4), Touch Probe (x2), Power Adapter, Protective Silicon Cover (x2), Spray Nozzle, O-Rings (x2), Disposable Sleeves (1 pack).

## Additional Accessories

Product	QTY per Pack	Item Code	Image
<b>MEASURING WIRE</b>	1 per Pack	EAPEXNMW	
<b>TOUCH PROBE</b>	2 per Pack	EAPEXTP	
<b>LIP HOOK</b>	2 per Pack	EAPEXLH	
<b>CONTRA ANGLE HANDPIECE CA 161</b>	1 per Pack	EAPEXCA	
<b>FILE CLIP</b>	4 per Pack	EAPEXFC	

**Warranty:** 2 years for handpiece motor, charging base, and power adapter. 1 year for contra angle attachment. 6 months for accessories.



# EDGEFLOW™

360° Rotating with excellent gutta-percha fluidity for a uniform, and smooth, bubble-free output ensuring seamless fillings.

## Benefits

- Excellent gutta-percha fluidity for a uniform, bubble-free output ensuring seamless fillings.
- Wireless and electric design for maximum comfort and ease of use.
- Dual battery system provides extended endurance, allowing for approximately 100 root canal fillings.
- Precise and intuitive display to better monitor the remaining gutta-percha.

### Specifications:

- Auto heating stop after timeout; Auto shutdown after timeout
- Charging base with battery
- Temperature monitoring design
- Disposable needles come pre-packed with gutta-percha
- Reusable injection needles available in various sizes and lengths
- Preset temperatures: 100°C, 120°C, 150°C, 180°C, 200°C
- Injection Speeds: Slow, Medium, Fast



## Accessories included in package

Item Code: **EEPTO**  
Main Unit: **EEPTORU**  
Charging Base: **EEPTOCH**

Main unit, Charging Stand, Power Adapter, Battery (x2), Heat Insulation Hood, Wrench of Needle, Gutta-Percha Ring, Injection Needle (x4) E23G-M (1) - E23G-L (1) - E25G-L (2).

## Accessories included in package

Product	QTY per Pack	Item Code	Image
<b>WRENCH</b>	1 per Pack	EEPTOW	
<b>LITHIUM-ION BATTERY</b>	2 per Pack	EEPTOB	
<b>THERMAL PROTECTOR CAP</b>	2 per Pack	EEPTOHIH3	

## Tips

Product	QTY per Pack	Item Code	Image
<b>*E20G-S INJECTION NEEDLE</b>	0 per Pack	EE20GSIN	
<b>*E20G-M INJECTION NEEDLE</b>	0 per Pack	EE20GMIN	
<b>*E20G-L INJECTION NEEDLE</b>	0 per Pack	EE20GLIN	
<b>E23G-M INJECTION NEEDLE</b>	1 per Pack	EE23GMIN	
<b>E23G-L INJECTION NEEDLE</b>	1 per Pack	EE23GLIN	
<b>*E25G-M INJECTION NEEDLE</b>	0 per Pack	EE25GMIN	
<b>E25G-L INJECTION NEEDLE</b>	2 per Pack	EE25GLIN	

Warranty: 1 year for device.

\*only available upon request



# EDGEPACK™

EDGEPACK™ Features high-precision temperature control with pre-set working temperatures, ensuring optimal compatibility with various gutta-percha types for precise and reliable obturation.

## Features



- Wireless Gutta-Percha Obturation Device for seamless, cordless operation.
- Ergonomic and ambidextrous design for enhanced comfort and usability.
- 0.35mm pre-bend tip for improved access and better penetration into narrow root canals.
- High-precision temperature control for consistent and reliable obturation.
- Efficient GP cutting system with high-conversion system.
- Compact size for easy cleaning and convenient storage.
- Visual and auditory prompts for enhanced safety during procedures.

### Specifications:


- Preset temperature options: 150°C, 180°C, 200°C, and 230°C.
- Advanced time-related features:
  - Ultra-fast warm-up time of 0.2 seconds.
  - Auto shut-off after 10 seconds for energy efficiency.
- Long-lasting battery life—provides up to 4 hours per charge.

## Accessories included in package









Item Code: **EEDPU**  
Main Unit: **EEDDURU**

Charging Stand, Power Adapter, WP4004 Tip, WP4504 Tip, WP5506 Tip, WP5508 Tip, Heat Insulation Sleeve of Tips (x4).

## Accessories included in package

Product	QTY per Pack	Item Code	Image
<b>CHARGING BASE</b>	1 per Pack	EEDPUCS	

## Tips

Product	QTY per Pack	Item Code	Image	Product	QTY per Pack	Item Code	Image
<b>WP4004 WORK TIP</b>	1 per Pack	EEDPUT4004		<b>EDGEPACK WP 5508 WORK TIP</b>	1 per Pack	EEDPUT5508	
<b>WP4504 WORK TIP</b>	1 per Pack	EEDPUT4504		<b>*EDGEPACK WP 5508L WORK TIP</b>	0 per Pack	EEDPUT5508L	
<b>WP5004 WORK TIP</b>	1 per Pack	EEDPUT5004		<b>*EDGEPACK WP 5510 WORK TIP</b>	0 per Pack	EEDPUT5510	
<b>*WP5506 WORK TIP</b>	1 per Pack	EEDPUT5506		<b>*EDGEPACK WP 6012 WORK TIP</b>	0 per Pack	EEDPUT6012	

Warranty: 1 year for device and tips.

\*only available upon request



# EDGEAPEX™

**EDGEAPEX™** Effectively locates the apical constriction, enabling precise determination of root canal working length for optimal treatment outcomes.

## Accessories included in package



Item Code: **EEALU**  
Main Unit: **EEALURU**  
Charging Base: **EEALUB**

### The EdgeApex is an apex locator and pulp tester.

It utilizes cutting-edge multi-frequency technology combined with a highly efficient digital signal processor, delivering optimal accuracy and information processing speed.

#### EDGEAPEX™ Apex Locator

- Charging Base (1)
- Measuring Wire (1)
- File Clip (4)
- Lip Hook (5)
- Touch Probe (2)
- Pulp Tester Probe (1)
- Tester (1)
- Lithium-ion Battery (1)
- Power Adapter (1)
- Endo Ruler (Available only upon request)
- Enlarged display of apical foramen (1)

## Additional Accessories

Product	QTY per Pack	Item Code	Image
<b>MEASURING WIRE</b>	1 per Pack	EEALUMW	
<b>TOUCH PROBE</b>	2 per Pack	EEALUMF	
<b>LIP HOOK</b>	5 per Pack	EEALULH5	
<b>ENDO RULER</b>	Available only upon request	EEALUER	
<b>FILE CLIP</b>	4 per Pack	EEALUFC5	
<b>TESTER</b>	2 per Pack	EEALUPTP	


**Warranty:** 1 year for device. 6 months for accessories.

# EDGEPOWER WORKFLOW

Explore our complete range of endodontic products on our website: [www.edgeendo.com](http://www.edgeendo.com)

## 1. ACCESS

- EDGEAPEX HP™
- EDGE R GLIDE
- EDGE APEX™ Locator




## 2. SHAPING

- EDGEAPEX HP™
- EDGEOne R




## 3. IRRIGATION

- Paper Points



## 4. OBTURATION

- EDGEBIOCERAMIC™ Sealer
- Gutta Percha Points
- EDGECore™
- EDGEFill™
- EDGEFlow™
- EDGEPack™



FOLLOW US ON  
SOCIAL MEDIA